

## Collecting EPA 549.2 Diquat

### Bottle Requirements

2- 250 ml Amber Glass  
preserved with Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>

### Sampling Location


The location code that is  
identified on the water system

### Sampling Plan.

### Sampling Instructions

1. Flush through the sampling point until main distribution water is coming through the sampling point (2-3 minutes).
2. Reduce the sample stream.
3. Fill and cap tightly and shake to ensure the preservative is dissolved.
4. Label the sample bottles and fill out the COC.
5. All samples must be placed on ice. If received next day the samples must be less than 6°C. If possible, deliver to the lab the same day that they are taken.

*Note: Holding Time is 7 days.*



**Accurate  
Environmental**

Engineering Services  
Environmental Training  
Laboratory Supplies  
Analytical Laboratories

**(800) 516 LABS (5227)**  
[www accuratelabs.com](http://www accuratelabs.com)

### 3 Locations in Oklahoma

#### **Accurate – Stillwater**

505 South Lowry Street Stillwater, OK 74074  
Fax (405) 372-5396

#### **Accurate – Tulsa**

3910 E. 51<sup>st</sup> Street Tulsa, OK 74135  
Fax (918) 663-6300

#### **Accurate – OKC**

12036 North Penn Ave. Okla. City, OK 73120  
Fax (405) 751-3108

*Need a quote or have questions?*

*Call our toll free number (800) 516-LABS  
or visit [www accuratelabs.com](http://www accuratelabs.com)*